

Colt Pursuit



Brand: Pride
Product Code: PCPURSUIT
Availability: 5
Weight: 97.50kg

Price: £3,889.00

Short Description

Features a powerful drivetrain and large 13" pneumatic wheels to handle rugged outdoor terrain. Luxury touches with a wraparound delta tiller.

Description

Luxury touches the Colt Pursuit scooter, like a wraparound delta tiller, a new LCD console and seat post suspension combine to give the Colt Pursuit a sleek, sporty style and a maximum level of comfort.

- Speed: 8mph
- Max Range: 21.8 miles
- Max Weight Capacity: 28.8 stone
- Ground Clearance: 9.69cm (3.81")
- Turning Radius: 170cm (67")
- Overall Length: 134.62cm (53")
- Overall Width: 63.65cm (25")
- Front and Rear Tyres: 10 cm x 33 cm (4" x 13") pneumatic
- Suspension: Front and rear
- Maximum obstacle climbing ability: 5cm (2")
- Total Weight without batteries: 97.5kg (215lbs)
- Heaviest Pieces when disassembled (rear section): 33kg (73.5lbs)
- Seat Weight: 21.77kg (48lbs)

- Seat Dimensions: 45cm wide x 45cm deep
- Drive system: rear-wheel drive, sealed transaxle, 24 volt DC motor
- Dual braking system: electronic, regenerative and electromechanical
- Battery Requirements: 2 x 12 Volt, deep cycle
- Battery Size: 50ah as standard (optional upgrade to 70ah)
- Battery Weight: 17kg or 38lbs for 50ah (24.27kg or 53.5lbs for 70ah)
- Battery Charger: 8amp off board

VAT RELIEF APPLIED
[Read more about VAT Exemption...](#)

Product Gallery



