

## Solo Bedstick Transfer



**Brand:** Aidapt  
**Product Code:** VY429/VY430  
**Availability:** 20  
**Weight:** 3.80kg

**Price: £44.55**

### Short Description

Solo Bedstick Transfer Aid assists transfer into and out of the bed. This effective transfer aid passes below the mattress and the looped handles extend up either side of the bed.

### Description

The Solo Bedstick transfer aid assists transfer into and out of the bed. The bedstick passes below the mattress and the looped handles extend up either side of the bed. The Bedstick transfer aid has a telescopic tube, so it can be safely adjusted to fit beds from 960 mm to 1580 mm. One side can be turned down for able bodied users sharing the bed, who do not require the transfer aid.

- High quality, durable finish
- Easy to stack
- Suitable for industrial cleaning
- Safety bar to prevent head entrapment
- 127 kg (20 stone) weight limit

	<b>VY429</b>	<b>VY430</b>
<b>Depth (mm)</b>	960-1580	960-1580
<b>Height (mm)</b>	420	420
<b>Width (mm)</b>	276	276

<b>Product Dimensions (mm)</b>	420x276x960	420x276x960
<b>Net weight (kg)</b>	3.8	3.8
<b>Colour</b>	White	Chrome
<b>Length (mm)</b>	960-1580	960-1580
<b>Maximum User Weight (kg)</b>	127	127
<b>Maximum User Weight (stones)</b>	20	20

**VAT RELIEF APPLIED**  
[Read more about VAT Exemption...](#)

## Product Gallery

