

Pride Go-Go Elite Traveller Plus



Brand: Pride
Product Code: PGG PLUS
Availability: 20
Weight: 39.20kg

Price: £1,789.00

Short Description

The first choice in travel mobility for many thousands of customers the world over. Truly transportable and passionate loved.

Description

Take the guess work out of travel mobility, with the Go Go Elite Traveller Plus.

With high performance operation and unbeatable ease of transport, the Go Go Traveller plus is sure to please.

Model Specification

Note: All specifications subject to change without notice.

- Weight Capacity: 23 Stone 3lbs
- Transport Weight: Total Weight (without batteries) 86.5 lb (39.2 Kg)
- Heaviest part (front section) 35 lb (16 Kg)
- Maximum Speed: 4 mph (6kph)
- Ground Clearance: 1.5" (3.8 cm)
- Braking System: Electronic, regenerative and electro - mechanical
- Turning Radius: 45.75" (116.2 cm)
- Standard Seating: Foldable, compact, vinyl-covered molded plastic

- Overall Length: 40.25" (102 cm)
- Battery Charger: Off-board, 2 amp
- Overall Width: 21.125" (54 cm)
- Battery Requirements: 2 x 12 volt - 18 amp (in one battery pack)
- Tyres: Front: 2" x 8" (5 cm x 20.32 cm) Rear: 3" x 9" (7.62 cm x 22.86 cm)
- Tip Wheels – Rear

[Read more about VAT Exemption...](#)

Product Gallery



