

Chelsfield High back chair



Brand: Aidapt
Product Code: VG818
Availability: 40
Weight: 0.00kg

Price: £351.60

Short Description

The Chelsfield Height Adjustable Chair is simple to maintain and clean.

Description

The Chelsfield Height Adjustable Chair is designed to offer comfort and longevity of service, as it is simple to maintain and clean for re issue.

The back of the Chelsfield Chair is ergonomically shaped for improved lumbar support and has a high back to offer head support and side head supports for comfort.

The frame is epoxy polymer coated and the arms are manufactured from polyurethane foam.

With fire retardant upholstery and height adjustable legs, this sturdy, wipe-clean chair is an attractive and long-lasting piece of furniture. The upholstery is now available in a range of colours to suit both modern and traditional interiors.

- Comfortable seating with lateral movement to suit the individual user
- High back for head support
- Fire retardant upholstery
- Specially designed legs to reduce tilt when transferring

- Easy clean tubular metal frame
- Simple height adjustment using 'e' clips
- 140 kg (22 stone) weight limit

Additional Specification

- Weight (kg): 14.55
- Seat Width (mm): 48.00
- Seat Depth (mm): 45.00
- Maximum User Weight (kg): 140.00
- Maximum User Weight (stones): 22.00

Product Gallery



