

## Card holder 2



**Brand:** NRS Healthcare  
**Product Code:** AA7311  
**Availability:** 30  
**Weight:** 0.00kg  
**Tax:** VAT 20% Applied

**Price: £9.95**

### Short Description

Allows those with reduced hand function to join in a game of cards.

### Description

Card holder 2 AA7311 is designed for people with limited use of their hands.

Cards can be easily added and removed. The front is extended to hold pencils or small coins.

### Specifications

#### Board:

- Height 35mm (13/8")
- Length 250mm (10")
- Width 50mm (2")