

Cobra



Brand: Drive Medical
Product Code: COBRA
Availability: 4
Weight: 98.00kg

Price: £5,460.00

Short Description

NEW stylish and distinctive 8mph Mobility Scooter is available now with a range of up to 32 miles.

Description

Drive Medical Cobra scooter is a NEW stylish and distinctive 8mph Mobility scooter available now with a range of up to 32 miles.

Features

- 4 wheel road legal mobility scooter
- 8mph (13kph) top speed
- Sliding and rotating seat
- Delta Bar
- High Level Indicators
- All round suspension
- 32 mile range
- 70ah batteries
- Stylish & distinctive design
- 29 stone user weight limit
- 12 months warranty

Specifications

- Length: 155cm (61")
- Width: 70.5cm (28")
- Max User Weight Capacity: 29 stone
- Range: up to 32 miles depending on user mass, terrain and battery options

[Read more about VAT Exemption...](#)

Product Gallery

