

Colt Executive



Brand: Pride
Product Code: COLTEX
Availability: 6
Weight: 117.80kg

Price: £4,409.00

Short Description

The Pride Colt Executive is the flagship of the Colt range of mobility scooters. A high performance 8mph road legal mobility scooter, the Colt effortlessly blends power and luxury.

Description

Colt Executive is the most stylish scooter, and comes with a wealth of high performance features as standard.

The Colt features all the high end features you could hope for: full suspension, sealed hydraulic brakes, extra large low-profile tires, advanced LED lighting, including directional signals, and quick adjust tiller.

The suspension features a mono-shock unit at the rear. This lets the Executive feel both comfortable and sporty, without the mushy feel of some large suspension scooters.

Plush captains style seat with a pillow top that is fully adjustable and rotates for ease of access.

An excellent choice for anyone looking for comfort and style backed up by Pride's renowned reliability.

- Wrap around delta tiller and LED lighting
- 25 - 30 miles range and user weight up to 28.5 stone
- Speed: 8mph
- Max Weight Capacity: 28.5 stone
- Ground Clearance: 13.33cm/5.25"
- Turning Radius: 210.cm/83"
- Overall Length: 142.8cm/56.25"
- Overall Width: 69.2cm/27.25"
- Front and Rear Tyres: 10 cm x 33 cm (4" x 13") pneumatic
- Suspension: Front and rear
- Total Weight without batteries: 117.8kg (259.lbs)
- Heaviest Pieces when disassembled (main frame): 93.21kg (205.5lbs)
- Battery Requirements: 2 x 12 Volt, deep cycle
- Battery Size: 100 AH (standard)
- Battery Weight: 27.44kg (60.5lbs) each
- Battery Charger: 8amp off board
- Insurance

VAT RELIEF
[Read more about VAT Exemption...](#)

Product Gallery



