

## Wheelchair crutch/walking stick holder and bag



**Brand:** Aidapt  
**Product Code:** VA132SS  
**Availability:** 30  
**Weight:** 0.75kg

**Price: £22.49**

### **Short Description**

Strong and spacious bag for a wheelchair or scooter. Features a zipped box lid and two side holders for your crutches or walking sticks.

### **Description**

The Wheelchair Crutch Bag has been specially designed to fit neatly onto the rear of most standard wheelchairs with pram handles, and many electric scooters too, and provides a generous, zipped storage space for coats or jackets, shopping, medication, snacks, etc. It also features a zipped rear compartment and two deep pockets on the front, one with a flap closed with 'hook and loop tape'; these are ideal for your umbrella, crutch, walking stick or other items. The Wheelchair Crutch Bag is a practical and affordable wheelchair accessory and enables wheelchair users and their carers to safely carry more, while keeping their hands free. It is water resistant and features a large zipped top closure for ease of use. There is a sturdy carry handle on the top of the bag, fixed loops on either side, which easily slip over the pram handles on the wheelchair, and adjustable straps with buckles for flexible fastening options. Raised rubber feet on the base keep the bottom of the bag dry when it's resting on the floor.

- Carry handle, fixed loops for handles and straps with buckles

- Generous storage space for coats, shopping, medication, snacks, etc
  - Large zipped top closure plus long pockets for crutches, walking sticks, etc.
  - Fits standard wheelchairs with pram handles
  - Enables wheelchair users and carers to safely carry more
  - Raised rubber feet keep the base of the bag dry
- 
- Depth (mm): 15
  - Height (mm): 415
  - Width (mm): 365
  - Product Dimensions (mm): 415x365x15
  - Net weight (kg): 0.54
  - Colour: Black
  - Maximum Load (kg): 12

[Read more about VAT Exemption](#)

## Product Gallery

