

Bed rope ladder



Brand: NRS Healthcare
Product Code: F30098
Availability: 15
Weight: 0.00kg

Price: £14.99

Short Description

Enables users with good upper limb strength to reposition themselves in bed.

Description

Bed rope ladder is a simple rope ladder attaches under the legs or castors at the end of a bed and allows the user to easily and simply pull themselves upright into a sitting position without assistance.

Features comfortable, wipe clean plastic rungs and synthetic rope. Suitable for all bed sizes.

Specifications

- Rung Length: 200mm (8")
- Diameter: 25mm (1")
- Cord Length: 3m (9ft 10")
- A managed 24 volt 300MA charger is provided
- Maximum user weight limit 140kg (22½st)

[Read more about VAT Exemption](#)



