

## Three tier fabric rise and recline chair



**Brand:** Drive Medical  
**Product Code:** CLR31  
**Availability:** 90  
**Weight:** 0.00kg

**Price: £659.00**

### **Short Description**

This Single Motor Rise and Recliner features a new compact modern design and a popular three tier waterfall back design.

### **Description**

The Three tier fabric rise and recline chair features a new compact modern design and a popular three tier waterfall back design. It is fully padded for increased comfort, with a three-tiered waterfall back.

- Stylish design
- Side storage pocket
- Ergonomically contoured seat & back
- Quiet and simple operation

### **Product Features:**

- Stylish design to suit any decor Full length chaise
- Side storage pocket
- Ergonomically contoured seat & back and padded upholstery for added comfort and support
- One motor activates a distinct lift, rest and recline action It will rise and tilt users to their feet

- Quiet motor provides a smooth and discreet operation Comes with a simple to operate two button handset
- Quality assured, meeting all relevant fire retardant standards
- Safety key is built into hand controller to prevent accidental movement and entrapment
- Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems

**Specifications:**

**Dimensions:** (W x D x H): 72cm (28") x 85cm (33") x 103cm (40.5")

**Seat:** (W x D x H): 53cm (21") x 54cm (21") x 51cm (20")

**Weight Capacity:** 140kg (22st)

[Read more about VAT Exemption](#)

**Product Gallery**



