

Comfort Knight Pillow Raiser



Brand: Aidapt
Product Code: VG814
Availability: 20
Weight: 13.66kg

Price: £705.00

Short Description

Comfort Knight Pillow Raiser enables one person to be safely raised and lowered in comfort, independently of any other occupant of the bed.

Description

The Comfort Knight Pillow Raiser bed aid enables one person to be safely raised and lowered in comfort, independently of any other occupant of the bed. This sturdy pillow raiser is powered by an extremely quiet actuator and is covered with a purpose made contoured mattress, finished in vinyl. Cotton covers are available. Suitable for use on the majority of beds. Not suitable for use with an adjustable bed.

- Ideal for people who cannot lie flat in bed
 - Finished in easy clean vinyl. Traversing safety straps as standard for increased stability
 - Contoured mattress with removable cover for simple cleaning
 - Easy touch hand control. Neat unobtrusive power box
 - 140 kg (22 stone) weight limit
 - Super quiet actuator with automatic safety spine
-
- Depth (mm): 150
 - Height (mm): 720

- Width (mm): 655
- Product Dimensions (mm): 720x655x150
- Net weight (kg): 13.66
- Angle of Raise (degrees): 8x68
- Colour: Beige
- Maximum User Weight (kg): 140
- Maximum User Weight (stones): 22
- Width (Canvas Size) (mm): 670

VAT RELIEF APPLIED
[Read more about VAT Exemption...](#)

Product Gallery

