

## Deluxe Multi Purpose Wheeled Table



**Brand:** Aidapt  
**Product Code:** VG832A  
**Availability:** 20  
**Weight:** 9.90kg

**Price: £103.96**

### Short Description

Deluxe Multi Purpose Wheeled Table is Ideal for those confined to bed and can be used around the home.

### Description

The Aidapt deluxe Multi Purpose Table is a convenient, mobile table for use around the home, particularly useful as an overbed table. This height-adjustable table is finished in an attractive beech melamine with genuine hardwood upstands and the top can also be set to different angles from flat to 45 deg. approximately, as required. In addition to the tilting table, the small side table provides excellent space for a drink or other smaller items.

- Tilting table with smaller side table for drinks etc.
- Ideal for use around the home
- Ideal for those confined to bed
- Maximum weight limit of 20kg.
- Height adjustable to suit most users.
  
- Depth (mm): 383
- Height (mm): 760-1170
- Width (mm): 808
- Product Dimensions (mm): 808x383x760

- Net weight (kg): 9.9
- Clearance required under bed (mm): 95
- Colour: White/Wood
- Folded Dimensions (mm): 885x140x420
- Footprint (WxD) Min/Max (mm): 600x400
- Height of Tube with Wheels: 65
- Overall Depth (mm): 808
- Overall Height (mm): 760-1170
- Overall Width (Folded) (mm): 140
- Overall Width (mm): 383
- Size: 808x383x760-1170

## Product Gallery

