

Solo Bed Transfer Aid



Brand: Aidapt
Product Code: VY428/ VY438
Availability: 15
Weight: 3.00kg

Price: £64.35

Short Description

Assists transfer into and out of one side of the bed, without affecting other occupants in the bed.

Description

The Solo Bed Transfer Aid assists transfer into and out of one side of the bed, without affecting other occupants. The hand grip sits above the bed and can be adjusted in height to suit, while the base frame is large enough to provide stability. The hand grip has been developed to MHRA guidelines. The Solo Bed Transfer Aid is available with or without a Strap and Buckle set for use on slatted beds.

- 7 height adjustment settings
- High quality easy to clean epoxy resin finish
- Large base for improved stability
- Height adjustable grip with safety rail

	VY428	VY438
Depth (mm)	680	680
Height (mm)	580-700	580-700
Width (mm)	600	600
Product Dimensions	580x600x680	580x600x680

(mm)		
Net weight (kg)	3	3
Colour	White	White
Configuration	Transfer Aid Only	Transfer Aid & Strap set
Handle Width (mm)	325	325
Maximum User Weight (kg)	127	127
Maximum User Weight (stones)	20	20
Overall Depth (mm)	680	680
Width of Base (mm)	600	600

VAT RELIEF APPLIED
[Read more about VAT Exemption](#)

Product Gallery

