

Adjustable Bed Cradle



Brand: Aidapt
Product Code: VY440
Availability: 20
Weight: 2.54kg

Price: £34.65

Short Description

Height Adjustable Bed Cradle is a Metal frame designed to keep bed covers away from the body.

Description

The Height Adjustable Bed Cradle from Aidapt is a metal frame designed to keep bed covers away from the body. Useful for those suffering from ulcers or skin irritation. The base is mounted between the mattress and the bed base and keeps bed clothes away from the feet and legs. Suitable for use at the side or end of the bed. Height is adjustable from 460mm to 670mm.

- Ideal for users with skin irritation or painful legs and feet
 - Durable construction and easy-clean finish
 - Can be fitted to the side or end of the bed as required
 - Height adjustable for a comfortable fit to suit each individual
-
- Depth (mm): 20
 - Height (mm): 355
 - Width (mm): 690
 - Product Dimensions (mm): 355x690x20
 - Net weight (kg): 2.54
 - Colour: White

- Overall Depth (mm): 360
- Width of Base (mm): 565

[Read more about VAT Exemption](#)