

## Triwalker with Seat



**Brand:** Drive Medical  
**Product Code:** TW017  
**Availability:** 30  
**Weight:** 0.00kg

**Price: £169.00**

### Short Description

The innovative and unique Triwalker with Seat provides a useful solution, enabling users to rest when needed.

### Description

With an attractive, metallic red steel frame, the Triwalker with Seat is lightweight and portable, offering stability and security to the user. Handle height can be easily adjusted and loop lock cable brakes can be easily applied.

- Unique seat design enables user to rest
- Loop lock cable brakes
- Lightweight frame weighing only 7kg (15.4lb)
- Handy side pockets for carrying essentials

### Product Features:

- Lightweight design at only 7kg (15.4lb)
- Unique seat design enables user to rest
- Front wheel swivels for increased manoeuvrability
- Easily folds for compact storage
- Easily adjustable handle height
- Max user weight - 114 kg (18 stone)

- 12 month warranty
- Includes bag, which folds with triwalker if stored away
- 7.5" wheels
- Suitable for indoor and outdoor use
- Easy to manoeuvre, durable steel frame

|                               |               |
|-------------------------------|---------------|
| <b>Specifications:</b> HEIGHT | 82-99cm       |
| WIDTH                         | 64cm          |
| DEPTH                         | 57cm          |
| HEIGHT MIN                    | 82cm          |
| HEIGHT MAX                    | 99cm          |
| WEIGHT                        | 7kg (15.4lbs) |
| SEAT WIDTH                    | 44cm          |
| SEAT DEPTH                    | 28cm          |
| SEAT HEIGHT MAX               | 53cm          |
| WEIGHT CAPACITY               | 114KG (18st)  |

[Read more about VAT Exemption...](#)

## Product Gallery





