

## Ambulance/transport 19 stone weight limit



**Brand:** Roma Medical  
**Product Code:** 1824  
**Availability:** 8  
**Weight:** 8.00kg

**Price: £262.00**

### Short Description

Well suited to moving patients in a confined space, the ambulance Transport or portorage wheelchair from Roma Medical is a very handy compact.

### Description

This ambulance chair folds down to an amazingly compact 9" deep when not in use and occupies just 48" in length making it very easy to collapse and store. The overall weight of the ambulance chair is just 7.5Kgs so can be easily lifted in and out of vehicles or storage.

### Features

- Easy-clean nylon upholstery
- 14" seat
- 4 Solid tyres
- Folds easily for storage / transportation
- Removable upholstery
- Lightweight for easy lifting

### Specifications

- Total weight: 7.5kg / 17lbs

- Maximum user weight: 19 Stone / 120Kgs
- Folded width 9" / 24cm
- Overall height: 40" / 101cm
- Overall width: 16" / 42cm
- Overall height: 39" / 99cm
- Wheelchair seat width: 14" / 35cm
- Wheelchair seat height: 20" / 51cm
- Wheelchair seat depth: 11" / 29cm

This chair has passed the crash test ISO7176-19:2001 for occupied transportation

\*A Crash Tested Wheelchair has a frame which is manufactured to comply with test ISO7176-19:2001 for occupied transportation

 [Read more about VAT Exemption...](#)