

Folding Commode



Brand: Drive Medical
Product Code: CO17
Availability: 20
Weight: 5.00kg

Price: £61.88

Short Description

Ideal folding commode for use when travelling or when space is limited. Comfortable armrests...

Description

Ideal folding commode for use when travelling or when space is limited.

Lightweight sturdy steel commode, folds easily for storage and transport. Comfortable armrests and a sort webbed backrest for added comfort.

Commode includes 5 litre potty with lid and handle that is mounted on a metal frame beneath the grey plastic seat.

The unique back of the fame forms a backrest when in use and a handle when folded.

Specifications

- Overall Width: 490mm (19")
- Height to moulded Seat: 495mm (19½")
- Total Height: 770mm (30")
- Seat Depth: 470mm (18½")
- Maximum User Weight Limit 115kg (18st)
- Product Weight 5kg (11lbs)

Product Gallery

