

Heavy duty 35 stone weight limit wheelchair



Brand: Roma Medical
Product Code: 1473
Availability: 6
Weight: 6.00kg

Price: £545.00

Short Description

Heavy Duty Self Propelled Wheelchair has a 22" seat width, Mag wheels, swing away footrests and padded armrests for added comfort. Maximum user weight 35 stone.

Description

The 1473 Heavy Duty Self Propelled Wheelchair comes with 24" pneumatic Mag rear wheels, height adjustable heel straps, parking brake, padded armrests, nylon upholstery, chromium plated steel and many more features as standard.

Features

- 60cm self-propelling mag puncture proof tyres
- 220cm front solid tyre swivel castors
- Double cross braces extra strength
- Dual axles for hemi-height convertability
- Detachable desk arms. Full arm optional
- Detachable swing away footrests with large footplate
- Elevating leg rest optional
- Nylon seat and back

Specifications

- Seat Depth 45cm
- Width Between Arms 56cm or 22"
- Depth 85cm
- Width 76cm
- Seat To Footrest 43-58cm
- Seat Height 53cm
- Height 94cm
- Weight 24kg
- Maximum User Weight 226kg or 35 stone

 [Read more about VAT Exemption.](#)