

Folding Seat Cane



Brand: Aidapt
Product Code: VP159A
Availability: 30
Weight: 0.00kg

Price: £23.09

Short Description

The Walking Cane with integrated Folding Seat is lightweight and completely portable.

Description

The Folding Seat Cane is a Walking Stick featuring a comfortable, curved foam handle, with a folding seat integrated into the frame. Manufactured from aluminium with a robust, plastic seat and non-slip feet, the tripod-style seat is suitable for users up to 16 stone (100 kg) in weight and weighs less than 1 kg, making it an ideal, portable solution.

- Foam handle
- Non-slip feet
- Lightweight, portable design
- Robust plastic seat
- Depth (mm): 60
- Height (mm): 850
- Width (mm): 280
- Product Dimensions (mm): folding:850x280x60, unfolding:
- Net weight (kg): 0.78
- Colour: Silver/grey
- Configuration: Cane with seat

- Design: SEAT CANE
- Height to Handgrip (mm): 850
- Maximum User Weight (kg): 100
- Maximum User Weight (stones): 15.7
- Number of Height Settings 1
- Overall Height (mm): 850
- Seat Height (mm): 510
- Size: 850x280x60

[Read more about VAT Exemption...](#)

Product Gallery

