

Heavy duty 32 stone weight limit wheelchair



Brand: Roma Medical
Product Code: 1485x
Availability: 6
Weight: 8.00kg

Price: £495.00

Short Description

The Roma 1485 heavy duty / bariatric wheelchair is a robust steel transit wheelchair. This wheelchair has a reliable and secure structure.

Description

The Roma 1485 heavy duty, bariatric wheelchair is a robust steel transit wheelchair. This wheelchair has a reliable and secure structure.

Features

- 30cm rear wheels
- 20cm x 50cm front puncture proof tyres swivel castors
- Double cross braces extra strength
- Dual axles for hemi-height convertibility
- Detachable arms
- Detachable swing away footrests with large footplate
- Elevating leg rest optional
- Nylon seat and back

Specifications

- Seat Depth 45cm

- Depth 71cm
- Width Between Arms 56cm
- Width 66cm
- Seat To Footrest 40-45cm
- Height 94cm
- Weight 20kg
- Seat Height 51cm
- Maximum User Weight 204kg

 [Read more about VAT Exemption...](#)