

Pvc Handledelbow Crutches



Brand: Aidapt
Product Code: VP147/8
Availability: 30
Weight: 0.00kg

Price: £39.99

Short Description

This lightweight aluminium elbow crutch has a wedge-shaped PVC handle for improved grip, along with a non-slip rubber foot.

Description

Pvc Handledelbow Crutches also features double height adjustment for comfort and flexibility. Priced per pair of elbow crutches, please ensure that you select the correct size from the available options when placing your order..

- Lightweight aluminium elbow crutch ideal for use during rehabilitation
- Double height adjustment with 14 plus 3 height settings
- Soft touch handgrip for greater comfort
- Maximum user weight limit 160kg (25 stone)
- Rubber foot to prevent slips

Specifications

Code	VP147	VP148
Colour	Silver/Grey	Silver/Grey
Depth (mm)	155	155
Height (mm)	775-1055	870-1155
Height to Handgrip (mm)	575-815	680-925

Inverter	No	No
Length (Collapsed) (mm)	775	870
Maximum User Weight (kg)	160	160
Maximum User Weight (stones)	25	25
Net weight (kg)	1.4	1.4
Number of Height Settings	17	17
Overall Height (mm)	1055	1155
Product Dimensions (mm)	n/a	680x100x155
Size	Medium	Adult (Large)
Width (mm)	100	100

 [Read more about VAT Exemption](#)