

Aktiv X3 Deluxe Lite Aluminium



Brand: Aktiv Wheelchairs
Product Code: X3
Availability: 8
Weight: 8.00kg

Price: £350.00

Short Description

Aktiv x3 Deluxe Lite Aluminium is a deluxe lightweight durable wheelchair for occasional to full-time use.

Description

This wheelchair has a lightweight aluminium folding frame and is easy to lift, store and transport. The armrests flip back and can be removed easily as standard and have comfortable polyurethane armrest pads. The backrest folds down and the leg rests and rear wheels are removable by easy quick release mechanisms; therefore the chair can be quickly and easily disassembled for easy storage in a car boot..

This chair seat height can be adjusted and can be easily and economically converted from a self-propel to an attendant transit wheeled model if required. The chair is fitted with anti-tip wheels for added safety and a stepper to help an attendant to ease the chair up over kerbs.

Features

- Folding Aluminum frame
- Fold down back
- Multiple axle for seat height + depth adjustments

- Flip back and detachable armrest with PU arm pad
- Plastic side panel clothes guard
- Swing away leg rests with black plastic foot plate (Angle adjustable) & heel loop
- Aluminum brake
- 7" grey PU castors
- Solid Tyres
- Anti-tipper
- Choice of press button or buckle style safety belt
- Black nylon seat
- Black nylon tension back (8 adjustable Velcro straps) with nylon cover

Optional

- Elevating leg rests

Specifications

- Seat Width: 41cm (16"), 46cm (18"), 51cm (20")
- Seat Depth: 43cm (17")
- Seat Height: Adjustable 41cm (16") - 51cm (20")
- Backrest Height: 42cm (16.5")
- Total Weight: 15.9kg
- Transport Weight: 10.8kg
- Max Weight: 130kg / 20.5st

This chair has passed the crash test ISO7176-19:2001 for occupied transportation

*A Crash Tested Wheelchair has a frame which is manufactured to comply with test ISO7176-19:2001 for occupied transportation



Product Gallery

