

Aktiv X6 Childrens Aluminium



Brand: Aktiv Wheelchairs
Product Code: X6
Availability: 6
Weight: 12.00kg

Price: £550.00

Short Description

Aktiv X6 Childrens Aluminium is a paediatric modular lightweight wheelchair suitable for occasional to full-time use and is multi-adjustable.

Description

The Aktiv X6 is a paediatric modular lightweight wheelchair suitable for occasional to full-time use and is multi-adjustable. This wheelchair has a lightweight Aluminium folding frame and is easy to lift, store and transport. The armrests flip back and can be removed easily as standard and have comfortable polyurethane armrest pads.

The leg rests and rear wheels are removable by easy quick release mechanisms; therefore the chair can be quickly and easily disassembled for easy storage in a car boot.

The chair can be easily and economically converted from a self-propel to an attendant transit wheeled model if required. The chair is fitted with anti-tip wheels for added safety and a stepper to help an attendant to ease the chair up over kerbs. Footplates are height and angle adjustable as standard and are supplied with a one piece footplate which simply slides onto the standard footplates to give a no-gaps platform for children to place their feet. Push handles are height adjustable as standard. The backrest is tension adjustable and the upholstery is

breathable and padded for extra comfort.

Features

- Folding aluminium frame
- Height adjustable push handles
- Flip back and detachable armrests with height adjustable PU arm pad
- Plastic side panel clothes guard
- Swing away leg rest with black plastic footplate, calf strap and one piece footplate cover
- Black nylon seat with breathable cover
- Tension adjustable backrest with breathable, wipeable cover
- Angle adjustable backrest
- 6" Grey PU front castors
- Available in Self Propel model (22" rear wheels)
- Solid Tyres
- Height adjustable anti-tippers
- Lap Belt
- 12" or 14" seat widths available
- Total weight - 14kg
- Transport weight 9.4kg

Dimensions

- Seat Width: 30cm (12"), 35cm (14")
- Seat Depth: 35cm (14")
- Seat Height: Adjustable 43cm (17") - 46cm (18")
- Backrest Height: 35cm (14")

Options

- Elevating leg rests

Colour available: Metallic Purple

This chair has passed the crash test ISO7176-19:2001 for occupied transportation

*A Crash Tested Wheelchair has a frame which is manufactured to comply with test ISO7176-19:2001 for occupied transportation

Product Gallery

