

Strood Height Adjustable Shower Stool



Brand: Aidapt
Product Code: VB635
Availability: 25
Weight: 2.80kg

Price: £54.45

Short Description

A high quality shower stool designed for heavy duty use; adjustable in height and featuring sturdy, integral side handles to aid transfer.

Description

The Strood Shower Stool has the added benefit of a clip-on/clip-off seat, making it simple to clean and reissue if required. The water outlets of the Strood fully conform to the MHRA guidelines, 1990 and 2004.

- Ideal for use in showering areas for people who struggle to stand for long periods
- Suitable for use in the majority of shower areas
- Height adjustable to suit most users
- Integral hand grips for extra security and a clip on/off seat for ease of cleaning

Specifications

- Lightweight and portable.
- Arm Height (mm): 560-690
- Colour: White
- Depth (mm): 420

- Footprint (WxD) Min/Max (mm): 515x420
- Height (mm): 560-690
- Height of Handle from Floor (mm): 560-690
- Maximum User Weight (kg): 159
- Maximum User Weight (stones): 25
- Net weight (kg): 2.8
- Overall Depth (mm): 420
- Overall Height (mm): 560-690
- Product Dimensions (mm): 560-690x515x420
- Seat Depth (mm): 275
- Seat Height (mm): 490-620
- Seat Width (mm): 370
- Width (mm): 515
- Width between Arms (mm): 465

 [Read more about VAT Exemption.](#)