

## Pride Apex Rapid



**Brand:** Pride  
**Product Code:** APEXRAP  
**Availability:** 60  
**Weight:** 39.60kg

**Price: £1,789.00**

### Short Description

Pride Apex Rapid is featuring all round suspension mobility scooter.

### Description

Pride Apex Rapid? with luxury and convenient features like a wraparound delta tiller with a high level charging port, LED kerb light and a storage basket together with a compact, easy to transport design.

Sleek styling complements the Apex Rapid's abundant features making it an exceptional travel mobility choice.

### Specifications:

|                                    |                                   |
|------------------------------------|-----------------------------------|
| Battery Weight (each)              | Weight: 18 AH: 13 kg (29 lbs.)    |
| Weight Capacity                    | 21 Stone 6 lbs (136 kg)           |
| Maximum Speed                      | 4 mph (6 km/h)                    |
| Ground Clearance                   | 3.8 cm (1.5")                     |
| Overall Length                     | 104.1 cm (41")                    |
| Overall Width                      | 52.1 cm (20.5")                   |
| Turning Radius                     | 116.2 cm (50.5")                  |
| Standard Controls                  | LED Display                       |
| Heaviest Piece (when disassembled) | 16.6 kg (36.6 lbs.) front section |

|                                  |                                                                                    |
|----------------------------------|------------------------------------------------------------------------------------|
| Total Weight (without batteries) | 39.6 kg (87.4 lbs)                                                                 |
| Suspension                       | Front and rear independent CTS Suspension                                          |
| Battery Requirements             | Two 12-volt, deep cycle (standard) Size: 2 x 18 AH                                 |
| Seating                          | Foldable, compact, vinyl-covered molded plastic                                    |
| Tyres                            | Front: 5 cm x 17.8 cm (2" x 7") solid<br>Rear: 6.35 cm x 20.3 cm (2.5" x 8") solid |

**VAT RELIEF**  
[Read more about VAT Exemption...](#)

## Product Gallery



