

Pride Colt Deluxe



Brand: Pride
Product Code: PC DELUXE
Availability: 19
Weight: 36.50kg

Price: £3,049.00

Short Description

6mph road/path friendly mobility scooter.

Description

The sleek, sporty Pride Colt Deluxe mobility scooter delivers smooth performance and all-new features.

With a full complement of advanced standard features including, rugged suspension, high-visibility automatic front and rear lighting system, feather-touch disassembly, LED kerb light, and wraparound, delta tiller, the Colt Deluxe electric disabled scooter is a greater value than ever before.

Front and rear suspension.

High-visibility automatic front and rear lighting system.

Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system.

Long-lasting, high-intensity LED headlight, tail light and indicators.

Pride's exclusive stylish, light-weight, non-marking, low profile wheels

Auto-connecting front-to-rear harness.

Delta tiller with wrap-around handles.

New back-lit battery meter for increased accuracy and visibility in any setting.

Includes pinchless hinges & eight- position seat rotation.

Easy-grip tiller adjustment knob.

Easy access tiller-mounted charger port for standard dual voltage off-board charger.

Rear view mirror standard

Specification

- Adjustable speed
- Front tyre size 9cm x 25cm
- Full suspension
- Ground clearance 3.5 cm
- Heaviest component 32
- Length 118.75 cm
- Maximum range 14 miles
- Maximum speed 6.25 mph
- Maximum user weight 159 kg
- Number of wheels 4
- Rear tyre size 10cm x 25cm
- Road legal
- Turning radius 137.16 cm
- Weight without batteries 36.5 kg
- Width 56.5 cm



Product Gallery

