

Drive Medical Envoy 4mph



Brand: Drive Medical
Product Code: DE4
Availability: 7
Weight: 94.00kg

Price: £2,210.00

Short Description

Drive Medical Envoy 4mph of the new Envoy range of Scooters are excellent additions to the Drive range, with the performance usually associated with much larger models.

Description

Drive Medical Envoy 4mph of the new Envoy range of Scooters are excellent additions to the Drive range, with the performance usually associated with much larger models.

These Scooters include a whole host of excellent up to date features:

- Max range of 30 mile *
- Top speed of 4mph *
- Powerful 350W motor
- All round suspension
- Stylish silver wheel hubs
- Delta bars as standard
- Black puncture proof non-marking tyres
- Battery - 50ah 12v x 2

* subject to conditions of use.

Specifications

- Length 121cm (48")
- Height 102cm (40")
- Width 60cm (24")
- Weight capacity 180kg (25st)
- Weight 94kg (207lb)

[Read more about VAT Exemption...](#)

Product Gallery

